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# ON THE RIGHT TRACK

**BEGO SEMADOS® RS/RSX-IMPLANTS**  
QUICK & FUNCTIONAL · UNIVERSAL & SUCCESSFUL  
RELIABLE & PROVEN · VARIED & SOPHISTICATED  
100% GERMAN QUALITY PRODUCT  
SIMPLE & INTUITIVE · QUICK & FUNCTIONAL  
UNIVERSAL & SUCCESSFUL · RELIABLE & PROVEN  
VARIED & SOPHISTICATED · SIMPLE & INTUITIVE

## The right product for many indications in implant dentistry

BEGO Semados® RS/RSX implants are individual and are available with machined shoulder (RS) or microstructured shoulder (RSX) and platform switching, depending on the clinician's preference and the requirements of the patient treatment. The range is supplemented by individual CAD/CAM-fabricated prosthetic components and surgical guides. The exceptionally high technical standards ensure the greatest satisfaction—for implant dentists and patients in equal measure.

### Your benefits

- Simple to use—insertion with self-centering
- Fast—just a few rotations needed when inserting
- Well-thought out—cover screw included
- For the most stringent quality requirements—100% German development and manufacturing
- Optional—Navigated/template-guided preparation of the BEGO Semados® RS/RSX-Line with diameters 3.75–4.5 can be carried out using the BEGO Guide system
- Established—internal taper connection with 45° medium taper angle and internal hex anti-rotation protection
- Integrated platform switching in the RS/RSX-Line supports the preservation of crestal bone

### Technical details

- Conical implant design with rounded apex
- Self-tapping thread design with optimal cutting angle
- Diameters: 3.0 / 3.75 / 4.1 / 4.5 / 5.5 mm
- Lengths: 7 / 8.5 / 10 / 11.5 / 13 / 15 mm
- Fabricated from grade 4 commercially pure titanium TiPure<sup>Plus</sup> surface
- Internal taper on the abutments ensures no gaps
- RS/RSX 3.0 is suitable for use in narrow anterior gaps to replace teeth 12, 22, 32–42
- Diameter of 5.5 mm in length 7 mm, for cases with a relatively low bone height


# Instructions

## BEGO Semados<sup>®</sup> RS/RSX-Implant

### 1 Pilot marker (initial drill)


Mark the bone and carry out initial drilling for directional orientation. Check with paralleling post.

### 2 Depth drill 2.5 (optional for RS/RSX 3.0)


Drill to the full implant length. Drill stops can be used. Check with paralleling post.

### 3 Final depth drill (in relation to the implant diameter)


Final enlargement of the cavity. Drill stops can be used.

### 3a Final depth drill for compact bone (alternative procedure to Pt. 3)


Final enlargement of the cavity with the DD D1 depth drills. Drill stops cannot be used.

#### 4 Countersink (optional)


Widen the cortical bone.

#### 5 Screw tapper (manual/handpiece)

Precut a thread in the compact bone.


#### 6 Implant insertion – manually or with a handpiece

Manual insertion

Insert the implant into the bone with the ratchet or the insertion tool used with a handpiece. (tightening torque  $\leq 50$  Ncm).


Insertion with handpiece

#### 7 Placing the cover screw

Pick up the cover screw with the hexagon screwdriver, insert into the implant and hand-tighten ( $\leq 10$  Ncm).

## Drill sequence

Implant	Pilot marker	Depth drill 2.5	Depth drill RS/R SX 3.0	Depth drill RS/R SX 3.75	Depth drill RS/R SX 4.1	Depth drill RS/R SX 4.5	Depth drill RS/R SX 5.5	Countersink*	Screw tapper
RS/R SX 3.0	x	(x)	x					(x)	(x)
RS/R SX 3.75	x	x		x				(x)	(x)
RS/R SX 4.1	x	x		x	x			(x)	(x)
RS/R SX 4.5	x	x		x	x	x		(x)	(x)
RS/R SX 5.5	x	x		x	x	x	x	(x)	(x)

x = required

(x) = optional (required for D1 bone quality)

\* = up to marking

## Rotational speed

Product	Rotational speed max. rpm
• Pilot marker / depth drill	→ 800
• Countersink	→ 800
• Screw tapper	→ 15
• Implant insertion	→ 15 (at ≤ 50 Ncm)

## Notes

Please read the instructions for use of the implant system (REF 84322) before using the system components.

**Caution:** The development of high temperatures, e.g., as a result of insufficient cooling and/or excessive pressing forces during preparation, must absolutely be avoided when using BEGO Semados® preparation instruments. In general, ensure that the BEGO Semados® preparation instruments are applied gently and carefully.

**Caution:** If very high tightening torques (> 50 Ncm) develop, the implant must be removed and set aside, while maintaining sterility, in the primary packaging because of risk of damage to the implant or the ratchet. The implant bed must be prepared again as necessary regarding the depth, drill hole diameter, bone thread or countersink for the head.

## BEGO Implant Systems GmbH & Co. KG

Wilhelm-Herbst-Str. 1  
 28359 Bremen, Germany  
 Phone +49 421 2028-246  
 Fax +49 421 2028-265  
 E-Mail [info@bego-implantology.com](mailto:info@bego-implantology.com)  
[www.bego.com](http://www.bego.com)